




APG as a reflection of vascular status





Dr. P. Kleikers
Mrs. N. Daimiwal



Principles of the APG technique


- APG is the second derivative of the fingertip photoplethysmogram (PPG) waveform (see also at 'APG technique').
 - Non-invasive
 - Pattern of 5 waves visible
- Used to measure vascular stiffness and aging
 - Arteriosclerosis: changes in the arterial wall lead to stiffer arteries with increasing age





APG and aging & vascular stiffness

- APG wave signals can be used to determine characteristics of vascular aging
- The pattern of 5 waves (a,b,c,d,e) shows characteristic changes with increasing age and increasing arterial stiffness¹.
- Ratios of the five waves change with increasing vessel stiffness
 - The b/a ratio increases with age, reflecting increased arterial stiffness.
 - The c/a, d/a and e/a ratio decrease with aging².
 - An aging index $((b-c-d-e)/a)$ was developed to determine the vascular age using the APG signal^{2,3}.





APG and cardiovascular diseases


- **APG ratios are correlated with**
 - Arterial distensibility⁴
 - Atherosclerotic changes (plaque formation) in the carotid artery^{4,5}
 - Blood pressure³
 - Gender³
- **Changes in the APG signal are also correlated with the risk of cardiovascular disease:**
 - Framingham risk score⁶
 - Changes in lipid profile and cholesterol, smoking and lack of regular exercise⁷
 - Hypertension and structural changes of the heart⁸





APG: other applications

- APG can also be used for other indications, for example:
 - Characterizing the response of the vessels to vasoactive agents^{2,9}
 - Characterisation of heart-rate-variability^{10,11}
 - Migraine¹²





APG versus Pulse Wave Velocity (PWV)

- **APG** is a relatively new technique in the cardiovascular field.
- The gold standard for measuring arterial stiffness is Pulse Wave Velocity
 - Technique depends also on the pulse wave travelling through the vessels.
 - **PWV**= the velocity of the pulse wave to travel a given distance between two sites of the arterial system¹³ (figure 1).





Figure 1 Schematic representation of Pulse Wave Velocity Measurement from carotid and femoral arteries ΔT =time delay between the foot waves; D =distance travelled by the wave.



APG versus Pulse Wave Velocity (PWV)

- **Correlation between APG and PWV**
 - The b/a and d/a ratio and the aging index of APG are correlated to PWV¹⁴.
 - In an hypertensive population, both PWV and APG could indicate vascular stiffness and aging¹⁵.
 - In a normotensive population, PWV and the APG indices were strongly correlated with age¹⁶.
 - In chronic kidney disease patients, APG correlated well with aortic calcification index as a measure of arteriosclerosis¹⁷.





APG: conclusion

The APG signal measured at the finger reflects the subject's overall vascular status and can give valuable information about the possible existence of cardiovascular diseases which warrant further investigation by a medical doctor. APG has been shown to correlate well with the gold standard (PWV) and is increasingly being used in the cardiovascular fields to easily and non-invasively measure vascular status.





For more in-depth discussion, please refer to the scientific publications.

1: Homma, S., Ito, S., Koto, T. & Ikegami, H. The relationship between accelerated plethysmography, blood pressure and arteriolar elasticity. *The Japanese Society of Physical Fitness and Sport Medicine* 41, 98-107, doi:10.7600/jspfsm1949.41.98 (1992).

2: Takazawa, K. et al. Assessment of vasoactive agents and vascular aging by the second derivative of photoplethysmogram waveform. *Hypertension* 32, 365-370 (1998).

3: Hashimoto, J. et al. Pulse wave velocity and the second derivative of the finger photoplethysmogram in treated hypertensive patients: their relationship and associating factors. *Journal of hypertension* 20, 2415-2422, doi:10.1097/01.hjh.0000042887.24999.7b (2002).

4. Imanaga, I., Hara, H., Koyanagi, S. & Tanaka, K. Correlation between wave components of the second derivative of plethysmogram and arterial distensibility. *Jpn Heart J* 39, 775-784 (1998)



For more in-depth discussion, please refer to the scientific publications.

5. Bortolotto, L. A., Blacher, J., Kondo, T., Takazawa, K. & Safar, M. E. Assessment of vascular aging and atherosclerosis in hypertensive subjects: second derivative of photoplethysmogram versus pulse wave velocity. *American journal of hypertension* 13, 165-171 (2000).

6. Otsuka, T., Kawada, T., Katsumata, M. & Ibuki, C. Utility of second derivative of the finger photoplethysmogram for the estimation of the risk of coronary heart disease in the general population. *Circulation journal : official journal of the Japanese Circulation Society* 70, 304-310 (2006).

7. Otsuka, T., Kawada, T., Katsumata, M., Ibuki, C. & Kusama, Y. Independent determinants of second derivative of the finger photoplethysmogram among various cardiovascular risk factors in middle-aged men. *Hypertension research : official journal of the Japanese Society of Hypertension* 30, 1211-1218, doi:10.1291/hypres.30.1211 (2007).



For more in-depth discussion, please refer to the scientific publications.

8. Iketani, T., Iketani, Y., Takazawa, K. & Yamashina, A. The influence of the peripheral reflection wave on left ventricular hypertrophy in patients with essential hypertension. *Hypertension research : official journal of the Japanese Society of Hypertension* 23, 451-458 (2000).

9. Millasseau, S. C., Kelly, R. P., Ritter, J. M. & Chowienczyk, P. J. The vascular impact of aging and vasoactive drugs: comparison of two digital volume pulse measurements. *American journal of hypertension* 16, 467-472 (2003).

10. Elgendi, M., Jonkman, M. & Boer, F. D. in 2010 *The 2nd International Conference on Computer and Automation Engineering (ICCAE)*. 514-517

11. Takada, M., Ebara, T. & Sakai, Y. The Acceleration Plethysmography System as a New Physiological Technology for Evaluating Autonomic Modulations. *General Health Examination* 35, 373-377, doi:10.7143/jhep.35.373 (2008)



For more in-depth discussion, please refer to the scientific publications.

12 Komatsu, K., Fukutake, T. & Hattori, T. Fingertip photoplethysmography and migraine. *Journal of the neurological sciences* 216, 17-21 (2003).

13 Safar, M. E., Henry, O. & Meaume, S. Aortic pulse wave velocity: an independent marker of cardiovascular risk. *The American journal of geriatric cardiology* 11, 295-298 (2002).

14 von Wowern, E., Ostling, G., Nilsson, P. M. & Olofsson, P. Digital Photoplethysmography for Assessment of Arterial Stiffness: Repeatability and Comparison with Applanation Tonometry. *PloS one* 10, e0135659, doi:10.1371/journal.pone.0135659 (2015).



For more in-depth discussion, please refer to the scientific publications.

15 Bortolotto, L. A., Blacher, J., Kondo, T., Takazawa, K. & Safar, M. E. Assessment of vascular aging and atherosclerosis in hypertensive subjects: second derivative of photoplethysmogram versus pulse wave velocity. *American journal of hypertension* 13, 165-171 (2000).

16 Hong, K. S., Park, K. T. & Ahn, J. M. Aging Index using Photoplethysmography for a Healthcare Device: Comparison with Brachial-Ankle Pulse Wave Velocity. *Healthcare informatics research* 21, 30-34, doi:10.4258/hir.2015.21.1.30 (2015).

17 Inuma, J. et al. Relationship between acceleration plethysmography and aortic calcification index in chronic kidney disease patients. *Hong Kong Journal of Nephrology* 14, 48-53, doi:https://doi.org/10.1016/j.hkjn.2012.09.003 (2012).



Contact

Head office:
Oxfordlaan 55
6229 EV Maastricht
The Netherlands
T: +31 (0)43 388 5841
E: info@vitak.com
W: www.vitak.com

Laboratories:
Oxfordlaan 70
6229 EV Maastricht
The Netherlands

Head office



Laboratories

